# Free-running and Q-switched Nd:YAG laser system, laser amplifier and second harmonic generation

### Goals

- Set-up and optimize Nd:YAG laser system in free-running red. Find the optimal output coupler reflectivity for maximum output energy from the set of three mirrors. Perform measurements and comparisons of energetic, temporal, and spatial laser beam characteristics.
- Determine the amplification coefficient G of an optical single-pass Nd:YAG amplifier.
- Set-up and run Nd:YAG laser system in passive Q-switch mode. Measure energetic and temporal laser beam characteristics.
- Set the KDP crystal to achieve the generation of the second harmonic frequency in the Q-switch mode. Observe the change in second harmonic radiation intensity (conversion efficiency) on the shade with changing the angle of rotation of the KDP crystal.

## **Experimental setup**



Figure 1: Simplified experimental setup.

## Instructions

## A) Laser characteristics in multi-mode free-running mode

- 1. Set up the center of laser crystal into the red laser diode beam.
- 2. Set up the laser resonator by beam coupling.
- **3.** Place the optical wedge behind the output mirror of the resonator, adjust the reflections of the laser diode radiation to detectors D1, D2.
- **4.** Run the laser (according to additional oral tutor instructions) and tune it on maximum output energy. Adjust the beam profile, recorded on a luminescent card or photosensitive paper, so that it is as close to a circular cross-section as possible.
- 5. Read the oscillator source voltage (the conversion table is on the last page) and compute the energy stored in the capacitors (C = 100  $\mu$ F) representing pumping energy E<sub>n</sub>.
- 6. Find the threshold pumping energy  $E_{th}$  for free-running mode. Measure the laser output energy E as a function of pumping energy  $E_p$ . Determine the maximum output energy  $E_{max}$  for the free-running mode.
- 7. Change the output coupler (R = 7 %, 48 %, and 74 % at 1064 nm) and repeat point 6.

- 8. In all other measurements, continue with the optimal output mirror (the highest  $E_{max}$  value).
- 9. Record the beam profile at maximum output energy  $E_{max}$  on black photosensitive paper. Determine its area S necessary to calculate the surface energy density  $W_{max} = E_{max}/S$ .
- **10.** Observe the generated radiation temporal evolution on oscilloscope and estimate generation time  $\tau_{FR}$  for various pumping energies:
  - i) just above E<sub>th.</sub>
  - ii) for middle of interval used E<sub>p</sub>.
  - **iii)** for maximal E<sub>p</sub>
  - Note: For free-running mode, the  $\tau_{FR}$  is the total duration of pulse envelope.
    - Reflectance of the wedge prism is R = 4%.
    - Transmittance of the ground-glass on the pyroelectric detector is T = 62%.
    - Sensitivity of the pyroelectric detector is 15.8 V/J.

#### **B)** Laser amplification

- 1. Setup laser amplifier active element to optical axis of justified oscillator.
- **2.** Setup up the wedge prism at the output of the Nd:YAG amplifier and displace energy detector to a new position (see Figure 1).
- **3.** Measure the output amplified energy  $E_{out}$  as a function of pumping energy  $E_p$  and compute an amplification coefficient  $G = I_{out}/I_{in} = E_{out}/E_{in}$ . Use the values of  $E_p$  for which you measured the energy  $E = E_{in}$  in the previous measurement (task A, point 6).

#### C) Laser characteristics in Q-switched mode

- **1.** Insert the Cr:YAG crystal for passive Q-switching into the resonator (using altitudinal feed) to have the red laser diode guide beam in the center and along the crystal axis.
- **2.** Transfer the wedge prism back to the original position beyond the output coupler and rearrange the energy detector to its previous location.
- **4.** Set the laser in Q-switched mode to generate one pulse per one pumping pulse. Observe and record one oscilloscope screenshot of the time evolution of the Q-switched pulse.
- 5. Determine the average length, energy and peak power, and their deviations of one Q-switched pulse (use at least 10 measurements). Don't forget to also record the beam profile to determine the surface energy density. Note: For the generation of Q-switched pulses, the pulse length  $\tau_Q$  means the full width at half maximum (FWHM) of the pulse amplitude.

#### D) Second harmonic generation

- 1. Continue with Q-switched mode and turn on the laser amplifier.
- 2. Adjust the KDP crystal to achieve second harmonic frequency generation.
- **3.** Observe the green radiation corresponding to the second harmonic frequency on the shade and try to change its intensity by changing the angle of rotation of the KDP crystal.

#### **Requested results**

A) Graph of the dependence of the maximum output energy  $E_{max}$  on the reflectivity of the output mirror  $R_{OC}$ . For all tested mirrors, graph and table of the dependence of the output energy E and efficiency  $\eta = E/E_p$  on the excitation energy  $E_p$ . For the optimal output mirror, graph and table of the dependence of the pulse length  $\tau_{FR}$ , and the average power  $P = E/\tau_{FR}$  on the

excitation energy  $E_p$ . Include the figures of the time evolution of radiation and the energy density at the maximum energy  $W_{max} = E_{max}/S$  (in J/cm<sup>2</sup>).

- B) Graph of the dependence of the amplification coefficient G on the excitation energy of the laser oscillator  $E_p$  for the optimal output mirror.
- C) Measurement results of the length  $\tau_Q$ , energy  $E_Q$ , peak power  $P_Q$  (all with deviation), surface energy density  $W_Q$  of Q-switched pulses and time evolution recording (from an oscilloscope) of a Q-switched pulse. Compare the surface energy density of the pulse in Q-switched mode and in free-running mode.

## References

Yariv, A., Quantum Electronics, chapters 7.3 – 7.4 rare-earth lasers MEOS - Q-switch theory (pages 9,10) <u>https://www.repairfaq.org/sam/MEOS/EXP0578.pdf</u>

## Appendix

| Set voltage [-] | Real voltage [V] | Set voltage [-] | Real voltage [V] |
|-----------------|------------------|-----------------|------------------|
| 200             | 369              | 480             | 498              |
| 210             | 370              | 490             | 506              |
| 220             | 371              | 500             | 515              |
| 230             | 373              | 510             | 524              |
| 240             | 374              | 520             | 534              |
| 250             | 376              | 530             | 543              |
| 260             | 378              | 540             | 553              |
| 270             | 381              | 550             | 563              |
| 280             | 384              | 560             | 574              |
| 290             | 387              | 570             | 585              |
| 300             | 390              | 580             | 596              |
| 310             | 394              | 590             | 607              |
| 320             | 398              | 600             | 619              |
| 330             | 402              | 610             | 631              |
| 340             | 406              | 620             | 643              |
| 350             | 411              | 630             | 655              |
| 360             | 416              | 640             | 668              |
| 370             | 421              | 650             | 681              |
| 380             | 427              | 660             | 694              |
| 390             | 433              | 670             | 708              |
| 400             | 439              | 680             | 722              |
| 410             | 445              | 690             | 736              |
| 420             | 452              | 700             | 750              |
| 430             | 459              | 710             | 765              |
| 440             | 466              | 720             | 780              |
| 450             | 473              | 730             | 795              |
| 460             | 481              | 740             | 810              |
| 470             | 489              | 750             | 826              |

**Table 1:** Conversion table between the numerical setting of the source voltage and the real value of the excitation voltage of the laser oscillator source. Source capacitor capacity:  $C = 100 \mu F$ .